


ACCURATE CELL COUNTING IS CRITICAL TO ENSURE EXPERIMENTAL SUCCESS AND REPRODUCIBILITY FOR CELL-BASED ASSAYS AND CELL CULTURE APPLICATIONS.

IMAGE-BASED COUNTERS

Image-based systems use bright field or fluorescent microscopy to capture an image of the cells. Some systems operate on a flow-based imaging methodology where cells are drawn into a capillary and the cells are imaged and counted as they pass a detector. Cell viability can be calculated using dye exclusion methods, such as Trypan Blue. Software analyzes the images based on system specific parameters such as cell diameter, brightness, and circularity to determine the number of cells and cell viability.



MASTERING CELL COUNTING

COULTER COUNTERS

Coulter counters measure changes in electrical resistance to determine the number, volume and size of the cells in the sample. Some Coulter counters offer the ability to distinguish live cells from dead cells and cell debris.



FLOW CYTOMETERS

Flow cytometers are not dedicated cell counters and report on relative values, such as the percent of cells in a given sample that have specific properties. The volume of sample counted needs to be determined to calculate the absolute cell count. To accomplish this, samples need to be spiked with fluorescent counting beads as a control.



CALCULATING THE NUMBER OF CELLS IS REQUIRED FOR MONITORING CELL PROLIFERATION & VIABILITY, OPTIMIZING CELL CULTURE CONDITIONS, AND PREPARING CELL-BASED ASSAYS.

CELL COUNTING PROCEDURES NEED TO BE:

- Accurate
- Reproducible
- Cell type independent, i.e., take into consideration cell property variables such as cell clumping, cell shape, and cell populations of varying size


HEMOCYTOMETER

DIMENSIONS	AREA	VOLUME AT 0.1 MM DEPTH
1 x 1 mm	1 mm ²	0.0001 mL
0.25 x 0.25 mm	0.0625 mm ²	0.00000625 mL
0.20 x 0.20 mm	0.04 mm ²	0.000004 mL
0.05 x 0.05 mm	0.0025 mm ²	0.00000025 mL

1. Apply the cell suspension so that the sample fills the chamber.
Tip: Dilute the cell suspension if necessary. Cells should be uniformly distributed without clumping or overlap on the grid.
2. Using a microscope, count the number of cells in 4 squares.
Tip: The lower the concentration of cells, the more squares should be counted to reduce statistical errors.
Tip: Only count cells that are within the square or on the bottom or left lines; do not count cells that touch the top or right lines of the square.
3. Calculate the density of cells in the suspension.

CELL ANALYSIS


CELL BRIGHTNESS



Cells & Clutter Counted

Only Cells Counted


CELL SHARPNESS



Sharp & Fuzzy Cells Counted

Fuzzy Cells Not Counted


CELL SPOT AREA



2% Area - Viable

10% Area - Dead

DECLUSTER DEGREE




None

Low

Medium


MINIMUM CIRCULARITY



Value of 1 Fewer Cells Counted

Value of .05 More Cells Counted

SPOT BRIGHTNESS



Default Value 75 - Classified Viable

Default Value 90 - Classified Dead

AUTOMATED CELL COUNTING

★★★★★★★★XX

VS.

MANUAL CELL COUNTING

★★XXXXXXX

Cell density (cells/mL) = $\frac{(\text{Average number of cells counted per square}) \times (\text{Dilution factor})}{\text{Volume of square (mL)}}$

Cell viability can be calculated by treating the cell suspension with cellular dyes, such as Trypan Blue. Live cells do not take up Trypan Blue, dead cells will stain blue.

Cell viability (%) = $\frac{\text{Total viable cells (unstained)}}{\text{Total cells (stained + unstained cells)}} \times 100$

AUTOMATED CELL COUNTERS OFFER AN EFFICIENT AND RELIABLE WAY TO QUICKLY COUNT AND ANALYZE A CELL POPULATION.

MANUAL CELL COUNTING OFFERS AN ACCESSIBLE WAY TO DETERMINE THE CONCENTRATION OF CELLS IN A LIQUID SAMPLE, REQUIRING JUST A LIGHT MICROSCOPE AND HEMOCYTOMETER.

★PROs: FAST • ACCURATE • REPRODUCIBLE • ELIMINATES SUBJECTIVITY • ANALYSIS & GRAPHICAL DISPLAY • DATA STORAGE & EXPORT CAPABILITIES • SUPPORTS 21 CFR PART 11 ✖CONs: INITIAL COST FOR EQUIPMENT

★PROs: INEXPENSIVE • VERSATILE ✖CONs: TIME-CONSUMING • LABOR-INTENSIVE • VARIABLE RESULTS • LIMITED FEATURES • POOR REPRODUCIBILITY • SUBJECTIVE

COUNT AND SIZE PARTICLES WITH HIGH RESOLUTION



Multisizer 4e COULTER COUNTER

Accurate counts for all types of particles — independent of shape, color, or transparency

- Precise and discrete size distributions in number, volume and surface area in one measurement (range: 0.2 μm to 1600 μm)
- Ultra-high resolution, multiple channel analysis and accuracy with the Digital Pulse Processor (DPP)
- Ability to detect small changes in particle/cell size
- Measurement is not affected by particle color, shape, composition or refractive index

www.cellcountsolutions.com

CHARACTERIZED
by ingenuity.



© 2015 Beckman Coulter, Inc. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

BECKMAN COULTER Presents

MASTERING CELL COUNTING

MASTERING CELL COUNTING

Beckman Coulter Life Sciences is dedicated to empowering discovery and scientific breakthroughs. The company's global leadership and world-class service and support deliver sophisticated instrument systems, reagents and services to life science researchers in academic and commercial laboratories, enabling new discoveries in biology-based research and development. A leader in centrifugation and flow cytometry, Beckman Coulter has long been an innovator in particle characterization and laboratory automation, and its products are used at the forefront of important areas of investigation, including genomics and proteomics. For more information, visit beckmancoulter.com. Follow Beckman Coulter Life Sciences on Twitter @BCILifeSciences; Facebook: BCILifeSciences; and LinkedIn.



**AUTOMATED
CELL COUNTING**
★★★★★★★x

HEAD
TO
HEAD

**MANUAL
CELL COUNTING**
★★xxxxxxx

CELL COUNTING MADE EASY



Vi-CELL XR Cell Viability Analyzer

Fast, Efficient, Accurate

- 12-position auto sampler, for walk-away sample analysis.
- No dilution required if your samples are within the 50×10^3 to 1×10^7 cells/mL range.
- Customizable cell types allow analysis of many cell species.

www.cellcountsolutions.com

CHARACTERIZED
by ingenuity.



Vi-CELL is For Laboratory Use Only. Not for use in diagnostic procedures.
© 2015 Beckman Coulter, Inc. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service names used herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

SPONSORED BY: 

CUSTOM PUBLISHING FROM: **TheScientist**
EXPLORING LIFE, INSPIRING INNOVATION